

PURPOSE

- To compare different type of Microcrystalline cellulose (MCC) as an excipient for Continuous Direct Compression system (CDC)
- To clarify the key parameter for CDC system
- ✓ MCCs with different particle sizes, bulk densities, and shapes were used (Table 1)
 - ★ Include trace impurity information (max or average of 6~10lots) as reference information since nitrosamine-related substances are important for pharmaceutical.


	KG-1000	KG-802	UF711	PH-102
Shapes	Fibrous	Fibrous	Amorphous (Porous)	Amorphous
Particle size D50 (μm)	50	50	50	90
Bulk density (g/mL)	0.12	0.21	0.22	0.30
Repose angle (°)	57	49	42	42
SEM Image (×500)				
★ Nitrite (μg/g) Maximum value	N.D.	N.D.	N.D.	0.012
★ Nitrate (μg/g) Average value	0.112	0.106	0.109	0.082

Table 1. Properties of MCCs used In this study.

EXPERIMENT METHODS AND RESULTS




Figure 1. Image of CDC system (CRA-RIS SYSTEM).

No.	1	2	3	4
Formulation	KG-1000	KG-802	UF-711	PH-102
Acetaminophen (APAP)	40	40	40	40
Spray-dried lactose (Lactose)	41	37	37	37
MCC: CEOLUSTM	15	20	20	20
Croscarmellose sodium (CCS)	2	1	1	1
Silicon dioxide (SiO2)	1	1	1	1
Magnesium stearate (Mg-St)	1	1	1	1
Bulk density (g/mL)	0.416	0.433	0.440	0.470
Repose angle (°)	46	43	42	42

Table 2. Formulation and properties of powder mix.

Tableting conditions
Tool: CRA-RIS SYSTEM (Kikusui Seisakusho)
45 punches, Turntable rotation speed: 41.2 rpm
Force feeder rotation speed: 45 rpm
CRATER screen diameter: φ3.0 mm (KG-1000), φ2.0 mm (KG-802, UF-711, PH-102)
Tablet size: 180 mg, φ8 mm - 12R
Run time: 60 min
Sampling time: 0, 2, 4, 6, 8, 10, 15, 30, 45, 60 min
Compression Force: (Main compression) 12 kN, (Pre-compression) 6 kN

Loss-in-Weight (LIW) feeders*	
POLARIS2 (20)	POLARIS5 (20)
Lactose 37-41 wt%	APAP 40 wt%
CCS 1 wt%	SiO2 1 wt%
POLARIS4 (20)	POLARIS3 (12)
MCC 15-20 wt%	Mg-St 1 wt%

* Due to the limited number of LIW feeders, lactose and CCS, APAP and SiO2 are premixed and set in the LIW feeder.

Continuous blending / tableting

No.		1	2	3	4
Items	Target	KG-1000	KG-802	UF-711	PH-102
Tablet Weight RSD	<2.0%	1.0%	1.0%	1.3%	1.0%
API content RSD	<2.0%	0.81%	1.36%	1.29%	1.87%
Tablet Hardness	>50N	100 N	91 N	93 N	81 N
Tablet Hardness RSD	-	7%	8%	7%	10%
Friability	<0.20%	0.15%	0.17%	0.17%	0.22%
Disintegrating Time	<30 min	1.9 min	2.1 min	2.0 min	1.5 min
Dissolution rate After 15 min	>85%	100%	100%	100%	99.7%

Figure 2. Experiment procedure and data summary

➤ KG-1000 showed the best performance among all MCC, especially API contents RSD was good.

DYNAMIC FLOWABILITY : BASIC FLOW ENERGY (BFE)




Figure 3. Illustration of measurement principle of powder rheometer FT-4.

➤ There were good correlation between API contents RSD and BFE.




Figure 4. API contents RSD and BFE data by MCC type.

FEED RATE STABILITY AND LONG RUN TEST




Figure 5. Feed rate of each CEOLUSTM grades during the run time.




Figure 6. Tablet hardness at each sampling time.